


For Converting

Alginates  
Nonwovens

Doyen's PCS system cuts and automatically places dressings into four side seal pouches.

### Benefits


- Speeds of up to 300 products per minute
- Toolless changeover via stepless electronic control of speed & cut-off length
- 2 unwinds for continuous process (splice on the fly) decrease product damage and waste
- Tensionless feed of delicate materials
- Extended infeed for re-feed of products with faulty pouches, web over option
- Fully validated, GAMP and GMP compliant



# PCS - Pad Cutting System



## TECHNICAL INFORMATION



Model	Pad Cutting System
Product length	50mm - 1000mm max.
Product width	50mm - 300mm
Product thickness	up to 4mm
Pouch width	90 - 350 mm
Speed	up to 300 cuts / min (cutter) up to 25 meters / min (4SS)
Reel specification	76mm core 600mm diameter max. (Pad material) 450mm diameter (Packaging material)

**Europe**  
Doyen Medipharm Ltd  
Gateway 11 Business Park, Copper Smith Way,  
Wymondham, Norfolk, NR18 0WY, UK  
Tel: (+44) 01953 600078 - Fax: (+44) 01953 605069  
sales@doyenmedipharm.co.uk

**North & South America**  
Doyen Medipharm LLC  
4030 South Pipkin Road, Suite 102,  
Lakeland, Florida 33811, USA  
Tel: (+1) 863-683-6335 - Fax: (+1) 863-683-6857  
salesUS@doyenmedipharm.com

**Asia**  
Doyen Medipharm Ltd  
5008 Jalan Lunas, First Floor, Batu 2 1/2  
Kulim, Kedah Darul Aman, Malaysia  
Tel: (+604) 495 1619 - Fax: (+604) 495 1697  
Surren@DoyenMY.com



[www.doyenmedipharm.com](http://www.doyenmedipharm.com)